# **Prizmatix**


# FC2-LED, FC3-LED and FC5-LED High Power Fiber Coupled Multi-Wavelength LED Light Sources

Ver. 06

#### Introduction

The new fiber-coupled high power UV, Violet, Blue, Green and Red LED Multi wavelength light source modules are effective replacements of lasers and lamps in many applications, such as spectroscopy, illumination and curing.

These 2, 3 and 5 LED modules provide high power CW or pulsed power at fiber output. The LED current controller supports CW operation mode with precise power control of each channel. The external TTL inputs enable external triggering of each LED. 3-LED and 5-LED are ideal for use with various fiber optic spectrometers in continuous, strobe or external triggering measurement mode.



Currently the 3-LED and 5-LED products can include any combination of following wavelengths: 365nm, 385nm, 390nm, 395nm, 400nm, 405nm, 410nm, 415nm, 420nm, 425nm, 430nm, 435nm, 440nm, 445nm, 455nm, 460nm, 465nm, 470nm, 475nm, 480nm, 495nm, 500nm, 505nm, 515nm, 520nm, 535nm, 595nm, 615nm, 630nm 640nm, 670nm, 740nm, 760nm, 810nm, 830nm, 940nm and White. Other wavelengths are available on request.

### **Features**

- High Power
- Reciprocal SMA fiber connection
- Precisely adjustable power by 10 turns potentiometer with dial
- Independently controlled power of each LED
- Speckle free
- Long life (no lamp replacement required)
- TTL external modulation input for each wavelength
- Analog modulation optional
- USB link for LED power computer control optional
- Replacement of multi wavelength lasers
- FC8-LED and FC11-LED with 8 and 11 LEDs also available
- Optional fiber bundles with homogenizer and collimator

Main Office	European Sales Office	North America Sales Office			
Phone: +972-72-2500097	Phone: +44(0)77-9172-9592	Phone: +1-(248)-436-8085			
Fax: +972-72-2500096	Fax: $+44(0)20-7681-2977$	Fax: +1-(248)-281-5236			
sales@prizmatix.com	sales.europe@prizmatix.com	sales.usa@prizmatix.com			
P.O.B. 244 Givat Shmuel 54101. Israel					

# **Prizmatix**

### **Specifications:**

#### Power output:

FC2-LED, FC3-LED and FC5-LED power output is limited by the reciprocal fiber characteristics. Larger core and higher NA will increase maximum power output, and vice versa. Please contact us for power data for specific wavelength and fiber.

The table below shows an example of the CW output power at three sample wavelengths for 1500µm core POF fiber:

Fiber brand	Fiber type*1	Pout @ 365nm Power Typ.*2	Pout @ 405 nm Power Typ.*2	Pout @ 445 nm Power Typ.*2
Mitsubishi Rayon SH6001 Super Eska	POF	100 mW	190 mW	180 mW

<sup>\*1:</sup> POF – Polymer Optical Fiber. NA=0.5, Core diameter=1500 $\mu$ m, Fiber length ~ 1m.

TTL input frequency: DC – 10 kHz (Higher freq. drivers available – please inquire)

**Analog input modulation frequency:** DC – 10 kHz (Optional)

Connector for TTL / Analog input: BNC Input power supply: 24 VDC, 1 A

Power adaptor input: 100-240 VAC, 1 A, 47-63 Hz

**Dimensions:** 

**FC2-LED:** 174mm x 130mm x 197mm (W x H x D) without extrusions. **FC3-LED:** 174mm x 130mm x 197mm (W x H x D) without extrusions. **FC5-LED:** 241mm x 130mm x 197mm (W x H x D) without extrusions.

Power adaptor: 60mm x 35mm x 100mm (W x H x L)



### **Optional Accessories:**

**Fiber patch cords:** Various fiber optic patch cords are available for use with the 3-LED or 5-LED products. Most popular patch cords are 1000 / 1500 / 2000 micron core diameter Polymer Optical Fibers (POF) terminated by optical SMA connectors on both sides. Prizmatix can provide Stainless Steel tube at one end of the patch cord instead of the SMA connector. This configuration is more convenient in some applications.

**Collimator:** The output from optical fiber is divergent according to fiber NA. In order to reduce the divergence angle a collimator module can be used. Prizmatix collimator was especially designed to fit thick core high NA Polymer Optical Fibers.

**Fiber Bundles:** To combine outputs of multiple LEDs a Y-shaped fiber bundle with two or more input branches can be used. Prizmatix can help to configure and build custom fiber bundles for specific applications. The output can be connected to homogenizer to mix all wavelengths into one homogenous beam.

**Special Spectroscopy Bundles:** For spectroscopy application special fiber optic bundles are very important. Special care must be addressed to eliminate possible cross link between the excitation and the collection fibers.

Main Office	<b>European Sales Office</b>	North America Sales Office			
Phone: +972-72-2500097	Phone: +44(0)77-9172-9592	Phone: +1-(248)-436-8085			
Fax: +972-72-2500096	Fax: +44(0)20-7681-2977	7 Fax: +1-(248)-281-5236			
sales@prizmatix.com	sales.europe@prizmatix.com	sales.usa@prizmatix.com			
P.O.B. 244 Givat Shmuel 54101, Israel					

<sup>\*2:</sup> Measurements were performed by Ophir Nova II power meter with PD300-UV head. The wavelength was set to 365 nm, 405 nm and 445 nm for each wavelength separately.