

FIBER COUPLED LED

Prizmatix


FC-BBW50 FC-BBW50-Plus



Fiber Coupled Broadband White LEDs

The FC-BBW50 and FC-BBW50-Plus feature wide visible emission spectrum. This broad spectrum is excellent for various reflection and transmission spectroscopy applications as well as color measurements. These self-contained light sources include all necessary driver electronics and thermal management. No external controller box is required just the included DC power adaptor.

- Ideal for spectroscopy applications & color measurements
- FC-BBW50 standard emission from ~ 410 to 670 nm
- FC-BBW50-Plus Extended emission from ~ 390 to 750 nm
- Self-contained – no controller needed
- Long Life LEDs
- Small footprint



Contact Prizmatix for information on High NA (0.63) Fiber Optics and Collimators ideal for LEDs.

Tel 248-436-8085
sales.usa@prizmatix.com
www.prizmatix.com

Prizmatix

LED Light Sources and Optical Systems